



PHOENIX PROGRAM

As Close As You Can Get To Bringing Discontinued Lenses Back To Life

Alden Optical's Phoenix Program offers eye care practitioners a simple solution to the problem of discontinued, terminally backordered, or capricious vial sphere lenses. Our unique blend of custom manufacturing, extraordinarily broad parameter ranges, and multiple lens material options means we're able to provide practitioners with Alden HP or Alden Classic lenses that are an extremely close substitute for many popular, but discontinued, sphere lenses. Contact us; chances are we can help.

Three Easy Steps to Transition Out Of Discontinued Lenses

1. Contact Alden Optical or an Alden Optical Authorized Distributor
2. Specify the habitual lens brand, specifications, and Rx
3. We will manufacture an Alden Optical lens with the same parameters in the material that most closely matches the habitual lens

Industry-Leading Service, Warranty, and Pricing

- We will ship the lens within 3 business days*
- We offer a full (100%) 90-day warranty on all single lenses

* Please allow an additional 2-3 days for tinted lenses

ALDENHP49

Hioxifilcon B, 49% water, 15Dk (ISO)

- Superior durability
- Good for dry eye lens wearers

ALDENHP59

Hioxifilcon A, 59% water, 24Dk (ISO)

- Superior end-of-day comfort
- Use when higher Dk is desirable

CLASSIC 38

Polymacon, 38% water, 9Dk (ISO)

CLASSIC 55

Methafilcon A, 55% water, 17Dk (ISO)





PHOENIX PROGRAM

BRAND CONVERSION TABLE

Manufacturer	Brand	Mfg Status	Alden Replacement	Note
Ciba Vision	Cibasoft SoftColors	discontinued	ALDEN Classic 38 Sphere	Add Alden enhancement tint
Ciba Vision	Cibasoft Visitint	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
Ciba Vision	CSI Daily Wear	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
Ciba Vision	CSI Toric	discontinued	ALDEN HP49 Toric	Specify back surface
Ciba Vision	Durasoft 2 Colors	discontinued	ALDEN Classic 38 Sphere	Add Alden enhancement tint
Ciba Vision	Durasoft 2 Colors for Light Eyes	discontinued	ALDEN Classic 38 Sphere	Add Alden enhancement tint
Ciba Vision	Durasoft 2 LiteTint	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
Ciba Vision	Durasoft 2 Optifit Toric Colors	discontinued	ALDEN HP49 Toric	Add Alden enhancement tint
Ciba Vision	Durasoft 2 Optifit Toric Colors for	discontinued	ALDEN HP49 Toric	Add Alden enhancement tint
Ciba Vision	Durasoft 2 Optifit Toric Handling Tint	discontinued	ALDEN HP49 Toric	Add Alden enhancement tint
Ciba Vision	D2 Special Eyes Prosthetic	discontinued	ALDEN PROSTHETIC	Alden can match most designs
Ciba Vision	Durasoft 3 LiteTint	discontinued	ALDEN Classic 55 Sphere	Alden can match all parameters
Ciba Vision	Durasoft 3 Optifit Toric	discontinued	ALDEN HP49 or Classic 55 Toric	Add Alden enhancement tint
Ciba Vision	D3 Special Eyes Prosthetic	discontinued	ALDEN PROSTHETIC	Alden can match most designs
Ciba Vision	Hydrocurve 3 Toric	discontinued	ALDEN Classic 55 Toric	Alden can match all parameters
Ciba Vision	Hydrocurve II	discontinued	ALDEN Classic 55 Sphere	Alden can match all parameters, DW only
Ciba Vision	Soft Mate B	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
Ciba Vision	Torisoft	discontinued	ALDEN HP49 Toric	Similar front surface toric design
Ciba Vision	Air Optix/02Optix Custom	discontinued	ALDEN HP59 Toric	Alden can match all parameters, DW only
CooperVision	Cooper Clear	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	Cooper Clear, CV43, & CooperHT	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	CooperClear FW & CV43 FW	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	Edge III	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
CooperVision	Expressions Accents	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	Frequency 38	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
CooperVision	Hydron Mini	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
CooperVision	Hydron Zero 4	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
CooperVision	Proclear Sphere (vial)	discontinued	ALDEN HP59 Sphere	Alden can match all parameters
CooperVision	Silver 07	discontinued	ALDEN Classic 38 Sphere	Alden can match all parameters
CooperVision	Vantage Accents	discontinued	ALDEN HP49 Sphere	Add Alden enhancement tint
CooperVision	Vantage Sphere	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	Vantage Sphere Thin	discontinued	ALDEN HP49 Sphere	Alden can match all parameters
CooperVision	Vantage Thin Accents	discontinued	ALDEN HP49 Sphere	Add Alden enhancement tint

Alden Optical does not offer a recommended alternative for Ciba Vision's D3 Colors, Colorblends, Complements, or Illusions brands.
Brands listed may be trademarks or registered trademarks of their respective companies.

RECOMMENDED TINT CONVERSIONS

Ciba Vision SoftColors	Alden Optical Enhancer	Durasoft 2 Colors/LE	Alden Optical Enhancer	Expressions Accents	Alden Optical Enhancer
Aqua	Azure #3	Sky Blue	Blue #2	Blue	Blue #3
Royal Blue	Blue #3	Violet Blue	Blue #3	Hazel	Walnut #1
Evergreen	Jade #3	Jade Green	Jade #3	Blue Topaz	Azure #3
Blue	Blue #2	Aqua	Azure #3	Aqua	Aqua #2
Green	Green #2			Gray	Gray #2
Amber	Walnut #1			Brown	Walnut #1
				Green	Green #3
				Jade Green	Jade #3



ALDEN
OPTICAL