

Tints Available

Alden HP Sphere & Toric, Alden Classic Sphere & Toric
 Prices represent tint fee only, and are single lens warranted price.

Standard Tints



- Visi-Tint** (edge to edge light blue)
- Tint Fee **call for pricing**



- Standard Tint** (11 mm tint diameter, transparent)
- Tint Fee **call for pricing**
 - Colors Aqua, Azure, Blue, Brown, Gray, Green, Jade, Walnut, Yellow
 - Shades Light (#1), Medium (#2), Dark (#3)



- Standard Tint with Clear Pupil**
- Tint Fee **call for pricing**
 - Pupil Diameter 2.0 mm to 8.0 mm in 0.5 mm steps, 4.0 mm used as standard

Prosthetic Tints



- Black Pupil** (Opaque)
- Tint Fee **call for pricing**
 - Pupil Diameter 2.0 mm to 13.0 mm in 0.5 mm steps



- Black Pupil with Standard Tint**
- Tint Fee **call for pricing**
 - Pupil Diameter 2.0 mm to 13.0 mm in 0.5 mm steps



- Black Annular with Clear Pupil**
- Tint Fee **call for pricing**
 - Annular Diameter 5.0 mm to 13.0mm in 0.5 mm steps
 - Pupil Diameter 2.0 mm to 8.0 mm in 0.5 mm steps



Tint Selection Guide



Alden Optical Laboratories
 13295 Broadway
 Alden, New York 14004

web www.AldenOptical.com
 email Info@AldenOptical.com
 tele 800.253.3669
 fax 800.899.5612

© July 2003
 Alden Optical Labs.