

A P P E A R I N G T H I S F A L L

BAUSCH+LOMB

ARISE™

ORTHO-K SYSTEM

DESIGN AN ORTHO-K LENS IN A FEW CLICKS
WITH CLOUD-BASED SOFTWARE THAT:



Starts smart

Identifies initial custom fit from corneal topography images.



Stays smart

Recommends refinements during postwear assessment.




Integrates easily

Instantly connects to select topographers and syncs existing patient data.

GET FIRST INSIDE LOOK

Crafts and assists with custom fits and refinements



1

CAPTURE CORNEA

Upload corneal image from your topographer.

2

GET CUSTOM DESIGN

Software identifies initial fit using thousands of image data points and patient inputs, such as eye shape and age. Adjust as you see fit.

Also recommends:

- o Landing zone—symmetric or toric
- o Optical zone—smaller or larger
- o Lens diameter—based on HVID

3

REFINE FIT

Upload new corneal image following overnight wear. If a refit is required, consider software-recommended refinements or use your expertise to adjust the fit.

ORTHO-K PRECISION THAT'S CONVENIENT

Auto-updates

Adds new features automatically—no manual downloads or updates.

Access records and support from anywhere

View patient data, history, and orders and chat with consultants.

Order online

Custom lens orders placed and processed digitally.

Visit [bausch.com/orthok](https://www.bausch.com/orthok) to be one of the first to know when it's available.

IMPORTANT SAFETY INFORMATION

The risk of microbial keratitis and other potentially vision-threatening adverse events has been shown to be greater among users of extended wear contact lenses than among users of daily wear contact lenses. If a patient experiences eye discomfort, excessive tearing, vision changes, or redness of the eye, the patient should be instructed to immediately remove the lenses and promptly contact his or her eye care practitioner.

See detailed Important Safety Information at [bausch.com/Eye-Care-Professional/Orthokeratology](https://www.bausch.com/Eye-Care-Professional/Orthokeratology)

BAUSCH + LOMB

See better. Live better.